



Correct



Not Correct

## Troubleshooting Guide

Problem: No/Few events being seen during acquisition

1) Is the green run light on?

**a. Yes**

1. Is your tube selected?
  - Arrow next to tube name should be highlighted green
2. Try a new experiment
  - Templates can unexpectedly become corrupt
3. Is the instrument connected properly?
  - Sometimes a "9" event threshold occurs
    - a. Turn off computer and instrument then turn on computer first and then instrument (and UV laser for Fortessa 1)
4. There might be bubble in flow cell
  - Prime 2x with nothing on sit (air bubbles should be coming out)
5. If none of the above fixes the issue, proceed to the step below.

**b. No**

1. Check/Refill tanks
  - Make sure o-ring is in place and all lines are connected correctly
  - Debubble sheath tank filters and side filter
  - Waste line may not be plugged in properly
2. Is your tube cracked?
  - Get new tube
  - Is your tube a 5 ml BD falcon **polystyrene**?
    - a. Change your tube to Catalog # 352058, 3502052, 352008, 352235
3. Are samples filtered and free of floaties?
  - Filter ALL samples prior to running
  - If running for longer than 1 hour, it is recommended to run a 3 minute cleaning of bleach and water.
4. Are your cells at the correct dilution?
  - 1-5 x 10<sup>6</sup> cells/mL is suggested
    - a. If not at correct dilution, dilute out sample with buffer or sheath
5. Instrument may be clogged.
  - Run 100% hot bleach for 3 minutes on HIGH followed by 5 minutes of H2O on HIGH



If this fails, please call 6-3573 or visit the lab, R5503. You can also text Alyssa @ 859-743-3323 or find any one of the RFCF Lab Members.