

## Reverse Transcription

4/27/06

1. Measure 1 ug for all RNA samples. Bring up to 8.5 ul with DEPC H<sub>2</sub>O.
2. Make a master mix of 0.5 ul of DNAase (1 u/ul) and 1 ul of 10x DNase Reaction buffer per sample.
3. Add 1.5 ul of Master mix for every sample.
4. Let samples incubate at room temperature for 15 minutes.
5. Add 1 ul of 25 mM EDTA to each sample.
6. Incubate samples at 65C for 10 minutes (PCR Machine).
7. Store at -80C or continue to next step.
8. Calculate master mix of 1 ul of Oligo (DT) and 1 ul dNTP mix per sample.
9. Add 2 ul of Master Mix to each sample.
10. Run program OligoDT (70C for 10 min and hold at 4C).
11. Store at -80C or continue to next step.
12. Calculate Master Mix at 4 ul of 5x First Strand Buffer, 1 ul DTTT, 0.25ul RT enzyme, 2.75 ul Sterile DI water per sample.
13. Add 8 ul of Master mix to each sample.
14. Run in program cDNA (42C for 1 hour, 72C for 5 min, hold at 4C).
15. Dilute samples with 100ul sterile DI water.
16. Store at -80 or -20C.