



## FEATURED ASSAYS

### IL-18

#### IL-18: What Is It?

IL-18 (also known as Interferon-gamma Inducing Factor or IGIF) is an important pro-inflammatory cytokine that is involved in a variety of immune processes associated with infection, inflammation, and autoimmunity. IL-18 is synthesized as a precursor molecule (pro-IL-18) in cells such as dendritic cells, macrophages, and intestinal epithelial cells. When activated, IL-1 $\beta$  converting enzyme (caspase-1) cleaves this precursor allowing it to be released from the cell. In combination with IL-12, IL-18 has been shown to act on T helper (Th1) cells to induce their production of Interferon-gamma (IFN-g). <sup>(1)</sup> IL-18 enhancement of Th1 cytokine (IL-2, GM-CSF, and IFN-g) production, FAS ligand expression, and IL-2R alpha chain expression has also been reported in Th1 cells. <sup>(2)</sup> IL-18 binding protein (IL-18BP) is a naturally occurring antagonist of IL-18. IL-18BP can downregulate TH1 responses by binding to IL-18, thus reducing the production of IFN-g.

#### IL-18: How Do We Measure It?

The Diagnostic Immunology Laboratory offers IL-18 testing. The IL-18 Immunoassay performed in our laboratory is a solid-phase ELISA (enzyme-linked immunosorbent assay) that employs the quantitative sandwich enzyme immunoassay technique. Specimen requirements are whole blood collected in an SST or red top tube. The tube may be spun once the clot has formed and preferably two serum aliquots of 0.3 ml placed at -20° C. When ready to ship, aliquots should be packaged on dry ice in a Styrofoam packing container for delivery the next day. Alternatively, we can accept unspun samples if they are received in our laboratory within 24 hours of collection and are kept at room temperature.

## IN THIS ISSUE:

### Featured Assays

CXCL9.....	1-3
IL-18 .....	3-4

Bulletin Board .....	4
----------------------	---


**Cincinnati Children's Hospital  
Medical Center  
3333 Burnet Avenue  
Cincinnati, Ohio 45229  
513-636-4200**

**Cancer & Blood Diseases Institute  
Clinical Laboratories**  
[www.cchmc.org/CBDILabs](http://www.cchmc.org/CBDILabs)  
Ph: 513-636-4685  
Fax: 513-636-3861

## IL-18: Why Measure It?


Measurement of IL-18 can aid clinicians in the diagnosis and monitoring of several immunological disorders. High levels of IL-18 are observed in patients with certain immune deficiencies and autoinflammatory diseases that involve activation or dysregulation of inflammasomes, such as XIAP deficiency and mutations of NLRC4.<sup>(3, 4, 5)</sup> These patients have persistently high levels of circulating total and also free IL-18 which contributes to disease pathogenesis (see Figure 1). Notably, a clinical trial is open to study the efficacy of a recombinant human IL-18 binding protein, Tadekinig alfa, in treating patients with XIAP and NLRC4 mutations.

(<https://clinicaltrials.gov/ct2/show/NCT03113760>).



**Figure 1.** Patient A presented to CCHMC with XIAP deficiency. IL-18 and CXCL9 (see discussion of CXCL9 on the next page) levels were used to monitor disease activity over the course of about 12 months. During this time the patient did experience disease flares indicated by the red arrows. The patient was treated and CXCL9 levels approached reference levels (dashed lines) during times of HLH remission. The IL-18, however, remained elevated even after recovery, indicating baseline ongoing inflammasome activity.

High IL-18 levels are also observed in patients with systemic-onset juvenile arthritis (sJIA) and may be used to help differentiate sJIA from other pyretic diseases. Studies have shown that IL-18 concentrations have a much higher predictive value than other serological biomarkers when trying to distinguish sJIA from diseases such as acute lymphoblastic leukemia, severe infections, Kawasaki disease and JIA, since the clinical symptoms of sJIA are often found in all these diseases.<sup>(6, 7)</sup> IL-18 can be used to help differentiate the diagnosis and prevent delay in treatment. Once the diagnosis has been made, it can be a beneficial marker to track disease activity and severity during management (see Figure 2).



**Figure 2.** Patient B presented to CCHMC with a diagnosis of sJIA and macrophage activation syndrome (MAS). IL-18 was used to monitor the patient during attempted disease control with biologics over the course of several months. These treatments were unsuccessful and plans were made to prepare the patient for an allogeneic hematopoietic stem cell transplant (HSCT). The HSCT was performed in December 2014 (indicated by the red arrow). Subsequently, IL-18 levels approached reference levels and remained there (dashed line).

High levels can additionally be observed in patients with active hemophagocytic lymphohistiocytosis (HLH), with levels of IL-18 usually returning to normal once HLH is in remission (except in patients with XIAP deficiency or NLRC4 mutations). Notably, the IL-18: CXCL9 ratio may help differentiate patients with sJIA and MAS (Macrophage Activation Syndrome) from patients with familial HLH.<sup>(8)</sup>

Lastly, IL-18 is also considered to be a secondary diagnostic criterion for autoimmune lymphoproliferative syndrome (ALPS).<sup>(9)</sup>

#### References:

1. Nakanishi K, Yoshimoto T, Tsutsui H, Okamura H. Interleukin-18 regulates both Th1 and Th2 responses. *Annu Rev Immunol.* 2001;19:423-474.
2. Dao T, Ohashi K, Kayano T, Kurimoto M, Okamura H. Interferon-gamma-inducing factor, a novel cytokine, enhances Fas ligand-mediated cytotoxicity of murine T helper 1 cells. *Cell Immunol.* 173(2):230-235.
3. Wada T, Kanegane H, Ohta K, Katoh F, Imamura T, Nakazawa Y, Miyashita R, Hara J, Hamamoto K, Yang X, Filipovich AH, Marsh RA, Yachie A. Sustained elevation of serum interleukin-18 and its association with hemophagocytic lymphohistiocytosis in XIAP deficiency. *Cytokine.* 2014 Jan; 65(1): p. 74-8. doi: 10.1016/j.cyto.2013.09.007. *Epub* 2013 Sep 29.
4. Canna SW, de Jesus AA, Gouni S, Brooks SR, Marrero B, Liu Y, DiMattia MA, Zaal KJ, Sanchez GA, Kim H, Chapelle D, Plass N, Huang Y, Villarino AV, Biancotto A, Fleisher TA, Duncan JA, O'Shea JJ, Benseler S, Grom A, Deng Z, Laxer RM, Goldbach-Mansky R. An activating NLRC4 inflammasome mutation causes autoinflammation with recurrent macrophage activation syndrome. *Nat Genet.* 2014 Oct; 46(10): p. 1140-6. doi: 10.1038/ng.3089. *Epub* 2014 Sep 14.
5. Romberg N, Al Moussawi K, Nelson-Williams C, Stiegler AL, Loring E, Choi M, Overton J, Meffre E, Khokha MK, Huttner AJ, West B, Podoltsev NA, Boggon TJ, Kazmierczak BI, Lifton RP. Mutation of NLRC4 causes a syndrome of enterocolitis and autoinflammation. *Nat Genet.* 2014 Oct; 46(10): p. 1135-1139. doi: 10.1038/ng.3066. *Epub* 2014 Sep 14.
6. Xia Y, Cui P, Li Q, Liang F, Li C, Yang J. Extremely elevated IL-18 levels may help distinguish systemic-onset juvenile idiopathic arthritis from other febrile diseases. *Braz J Med Biol Res.* 2017 Feb 16; 50(2):e5958. doi: 10.1590/1414-431X20165958.
7. Takahara T, Shimizu M, Nakagishi Y, Kinjo N, Yachie A. Serum IL-18 as a potential specific marker for differentiating systemic juvenile idiopathic arthritis from incomplete Kawasaki disease. *Rheumatol Int.* 2015 Jan; 35(1): p. 81-4. doi: 10.1007/s00296-014-3059-2. *Epub* 2014 Jun 12.
8. Weiss Eric S, Girard-Guyonvarc'h Charlotte, Holzinger Dirk, de Jesus Adriana A, Tariq Zeshan, Picarsic Jennifer, Schiffrin Eduardo J, Foell Dirk, Grom Alexei A, Ammann Sandra, Ehl Stephan, Hoshino Tomoaki, Goldbach-Mansky Raphaela, Gabay Cem, Canna Scott W. Interleukin-18 diagnostically distinguishes and pathogenically promotes human and murine macrophage activation syndrome. *Blood* 2018 131: p. 1442-1455; doi: <https://doi.org/10.1182/Blood-2017-12-820852>.
9. Oliveira, JB, Bleesing JJ, Dianzani U, Fleisher TA, Jaffe ES, Lenardo MJ, Rieux-Laucat F, Siegel RM, Su HC, Teachey DT, and Rao VK. Revised diagnostic criteria and classification for the autoimmune lymphoproliferative syndrome (ALPS): report from the 2009 NIH International Workshop. *Blood.* 2010. 116(14): p. e35-40.

## CXCL9

### Measurement of CXCL9: What Does It Mean?

Chemokines are types of cytokines that cause various cells to move (chemokinesis), usually towards inflammation. There are numerous chemokines that act in the immune system, with a confusing array of receptors. CXCL9 is a chemokine, previously known as 'monokine induced by interferon gamma' (MIG), which binds the chemokine receptor CXCR3 along with two other chemokines, CXCL10 and CXCL11.<sup>1,2</sup> These 3 chemokines have a number of roles, but their main function is to attract T cells into inflamed tissues. As its prior name suggests, CXCL9 is secreted by monocytes and macrophages in response to another cytokine, interferon gamma (IFN- $\gamma$ ).

The CXCL9 Immunoassay performed in our laboratory is a solid-phase ELISA that employs the quantitative sandwich enzyme immunoassay technique. Specimen requirements are whole blood collected in an EDTA tube that has been spun within 8 hours of collection. Divide all available plasma into 0.5 ml aliquots. Three aliquots of 0.5 ml of plasma is preferred, one is acceptable. Freeze immediately at -20° C or below and ship to our laboratory on dry ice.

### CXCL9: How Do We Measure It?

CXCL9 is unique because its release is triggered almost entirely by IFN- $\gamma$  (IFN- $\alpha$  and TNF- $\alpha$  have a slight effect on it in some contexts).<sup>1,2</sup> Because of this specificity, it is an excellent biomarker for the action of IFN- $\gamma$ . While IFN- $\gamma$  can be measured directly in plasma, measurement of CXCL9 in experimental contexts appears to be a much more sensitive measure of IFN- $\gamma$  activity. One possible explanation of this is the fact that IFN- $\gamma$  is secreted by lymphocytes in a very directional fashion (a 'private message' for the interacting cell), while CXCL9 is secreted by mononuclear cells in all directions (a 'public broadcast' calling other cells towards inflammation).<sup>3</sup>

## What is known about CXCL9 in disease?

CXCL9 may be helpful when diagnosing inflammatory conditions where IFN- $\gamma$  is important. Specifically, IFN- $\gamma$  appears to play a significant role in hemophagocytic lymphohistiocytosis (HLH) and macrophage activation syndrome (MAS). Elevated CXCL9 is seen in experimental models of these diseases.<sup>4 5</sup> In patients, it is reported to be elevated in lymphoma-associated HLH<sup>6</sup> and MAS<sup>7</sup>, while our unpublished observations indicate that it is elevated in essentially all patients with HLH. Though these observations suggest that CXCL9 is very sensitive for detecting HLH/ MAS, it is not likely to be entirely specific for these conditions, as IFN- $\gamma$  is involved in many immune responses.

CXCL9 may also be useful for assessing response to IFN- $\gamma$  blocking therapies. An investigational agent which blocks IFN- $\gamma$  (emapalumab) is currently being studied as treatment for HLH and MAS, and JAK inhibitors (blocking IFN- $\gamma$  signaling) have been suggested as experimental therapy for these disorders. In both cases CXCL9 may be a useful marker of IFN- $\gamma$  blockade in patients and may prove to be a valuable part of routine therapeutic monitoring.

### References:

1. Groom JR, Luster AD. CXCR3 ligands: redundant, collaborative and antagonistic functions. *Immunol Cell Biol.* 2011;89(2):207-215.
2. Corbera-Bellalta M, Planas-Rigol E, Lozano E, et al. Blocking interferon gamma reduces expression of chemokines CXCL9, CXCL10 and CXCL11 and decreases macrophage infiltration in ex vivo cultured arteries from patients with giant cell arteritis. *Ann Rheum Dis.* 2016;75(6):1177-1186.
3. Huse M, Quann EJ, Davis MM. Shouts, whispers and the kiss of death: directional secretion in T cells. *Nat Immunol.* 2008;9(10):1105-1111.
4. Buatois V, Chatel L, Cons L, et al. Use of a mouse model to identify a blood biomarker for IFN $\gamma$  activity in pediatric secondary hemophagocytic lymphohistiocytosis. *Transl Res.* 2017;180:37-52 e32.
5. Prencipe G, Caiello I, Pascarella A, et al. Neutralization of IFN-gamma reverts clinical and laboratory features in a mouse model of macrophage activation syndrome. *J Allergy Clin Immunol.* 2018;141(4):1439-1449.
6. Maruoka H, Inoue D, Takiuchi Y, et al. IP-10/CXCL10 and MIG/CXCL9 as novel markers for the diagnosis of lymphoma-associated hemophagocytic syndrome. *Ann Hematol.* 2014;93(3):393-401.
7. Bracaglia C, de Graaf K, Pires Marafon D, et al. Elevated circulating levels of interferon-gamma and interferon-gamma-induced chemokines characterise patients with macrophage activation syndrome complicating systemic juvenile idiopathic arthritis. *Ann Rheum Dis.* 2017;76(1):166-172.

## Acknowledgments

The authors would like to thank **Michael Jordan MD**, Division of Immunobiology, CCHMC for his contribution of the CXCL9 portion of this Newsletter.

The Newsletter content writer is Sabina Sylvest  
The Newsletter is edited by Rebecca Marsh, MD

## BULLETIN BOARD

Our laboratory recently completed an 8-week course using hands-on exercises, team-based projects and 1-on-1 coaching to deepen each participant's ability to solve problems and improve operational performance. The course was based on the lean tools & methods used by Toyota and top performing healthcare organizations to identify non-value added work and lead improvements in safety, quality, productivity, lead time, and cost effectiveness. Cincinnati Children's Hospital is committed to empowering our frontline staff to drive change to eliminate non-value added work and improve our customer (you and your patients) satisfaction!

In our project, we tackled our specimen intake process. Receiving over 100 tests a day from our external clients is no small feat! By reviewing this process using our Lean methodology, we were able to reduce our specimen processing time by 37% without affecting quality. This allows us to get back to being available for our customers sooner, expedite initiation of sample testing, and reduce over-processing. We want you to know that we are committed to changing the outcome and improving our services. This was a great introductory project that allowed us to become familiar with these tools that we can then carry forward to other improvements within the laboratory.