



## Fall Meeting

Please stop by and say “hello” in December at this year’s American Society of Hematology (ASH) Annual Meeting in New Orleans. We will be there as the Cincinnati Children’s Clinical Laboratories.

### THIS ISSUE’S FOCUS – Th17 Assay

Hyperimmunoglobulin E syndrome (HIES) is a rare immune deficiency typically characterized by dermatitis, recurrent infections, abnormalities of bone and connective tissue and elevated serum IgE.<sup>1</sup> The most common genetic mutations reported in HIES patients are in either STAT3 or DOCK8 which result in interference with the differentiation pathway of Th17 cells, causing low or non-existent cell numbers.<sup>2</sup>

Th17 is the name of interleukin (IL)-17 producing CD4+ T cells that are a subpopulation of T helper lymphocytes involved in immune responses to fungal and extracellular bacterial pathogens. We have developed a flow cytometric assay that detects the peripheral blood levels of Th17 cells in order to screen patients for Hyper IgE Syndrome.

In our assay, IL-17 cytokine production is stimulated by activating whole blood in the presence of Brefeldin-A to block release of intracellular cytokines, trapping cytokines produced during activation inside the cell. The activated cells are stained with surface monoclonal antibodies for phenotyping, fixed, permeabilized and stained with antibody against IL-17. The cells are then analyzed by flow cytometry. See Figure 1 and 2 for examples of the analysis on a healthy adult control as well as on a patient with HIES.

### IN THIS ISSUE:

|                               |            |
|-------------------------------|------------|
| <b>Th-17: New Assay</b> ..... | <b>1-4</b> |
| <b>Bulletin Board</b> .....   | <b>5</b>   |
| <b>Current Menu</b> .....     | <b>6</b>   |

Cincinnati Children’s Hospital  
Medical Center  
Cancer & Blood Diseases Institute  
3333 Burnet Avenue  
Cincinnati, Ohio 45229  
513-636-4685

#### Director

Alexandra H. Filipovich, MD

#### Associate Director

Jack J.H. Blessing, MD, PhD


#### Director of Test Development

Rebecca A. Marsh, MD

#### Laboratory Supervisor

Darryl Hake

### Adult Control



### HIES Patient




Figure 1. Adult control showing normal expression of Th17 cells. Example analysis of the Th17 assay. We first classify cells by surface phenotyping of CD4+CCR6+CD45RA- expression. We then look at the isotype staining on activated cells as well as the IL17+ staining on non-activated cells to set our IL17+ gate. 38.8% of CD4+ cells are CCR6+CD45RA-. 3.5% of CCR6+CD45RA- are IL17+.

Figure 2. HIES patient showing decreased to absent IL-17 production, reflecting a lack of TH17 cells. Only 0.5% of the CD4+CCR6+CD45RA- cells are IL17+.

Regardless of the % of CD4+ T cells which were CD45RA- CCR6+, the % of Th17 cells (that produce IL-17) was fairly constant among healthy adult controls and ranged from 0.31-1.8% of the total CD4+ T cell population (Figure 3 and Table 1).

To establish normal pediatric ranges, peripheral blood from healthy children was used to determine pediatric normal ranges of T cell populations which are CD4+ CD45RA- and CCR6+ (Table 1). This information was then used in conjunction with the adult normal range for Th17/IL17A+ cells to screen for Hyper IgE Syndrome in pediatric patients.

A limitation of the assay is our ability to run samples at various timepoints post-collection. Comparisons revealed that surface marker expression does not change following activation of samples, but that surface marker/T cell phenotype differs in samples tested immediately after blood draw versus samples which were 24 hours old. Therefore, all testing must be performed on samples rested at room temperature for 24 hours following blood draw in order to have a viable reference range.

Table 1. Shows normal ranges established for each result component as well as a result from a confirmed HIES patient for comparison.

|                    | CCR6+CD45RA-<br>(of CD4+) | IL17+<br>(of CD4+CCR6+CD45RA-) | Th17<br>(of CD4+) |
|--------------------|---------------------------|--------------------------------|-------------------|
| Adults             | 15.6 – 47.0 %             | 1.2 – 4.6 %                    | 0.31 – 1.80 %     |
| <2 years           | 2.4 – 6.1 %               |                                |                   |
| 2-6 years          | 4.2 – 12.2 %              |                                |                   |
| 6-18 years         | 10.7 – 26.9 %             |                                |                   |
| HIES patient, 3 yo | 6.8%                      | 0.5%                           | 0.04%             |

### Patient Vignette:

A three-year old male was referred to the clinic for evaluation of recurrent infections and concern for Hyper IgE syndrome. His IgE was 628 (2-199 IU/mL ref range) with a history of rashes, presumed eczema, that did not respond to typical treatments, as well as respiratory infections. The Th17 assay performed in our laboratory revealed as absence of IL-17+ CD4+ CCR6+ CD45RA- cells, indicating a lack of TH17 cells (Figure 2 and Table 1). STAT3 sequencing was performed which showed a novel missense mutation (heterozygous A>G nucleotide substitution) and confirmed the diagnosis of HIES. Treatment continues to be prophylactic oral antibiotics as well as topical antibacterials and topical steroids.

### References

1. Paulson ML, Freeman AF, Holland SM. Hyper IgE syndrome: an update on clinical aspects and the role of signal transducer and activator of transcription 3. *Curr Opin Allergy Clin Immunol* 2008; 8:527.
2. Holland SM, DeLeo FR, Elloumi HZ, et al. STAT3 mutations in the hyper-IgE syndrome. *N Engl J Med* 2007; 357:1608.

Figure 3. Results from 40 healthy adult donors.





Figure 4. Surface Marker Expression in healthy pediatric donors



The Diagnostic Immunology Laboratories, consisting of the Clinical Immunology Laboratory and the Research Immunology Laboratory, are committed to providing the highest quality, comprehensive clinical testing available to aid in the detection, diagnosis and treatment of pediatric immunologic, as well as oncologic and hematologic, disorders. We're committed to applying scientific advances to promote efficiency, enhance patient care and improve clinical utility.

The clinical diagnostic laboratories are in compliance with all major regulatory agencies including CLIA (Clinical Laboratory Improvement Amendments), CAP (College of American Pathologists), HCFA (Health Care Financing Administration), HIPAA (Health Insurance Portability and Accountability Act) and JCAHO (Joint Commission on Accreditation of Healthcare Organizations).

The current menu of immunologic assays and information regarding shipping instructions is published on the last page of this Newsletter. The accompanying Test Requisition Form can be obtained through our website. Previous editions of the Newsletter can also be found at this website: [www.cchmc.org/DIL](http://www.cchmc.org/DIL)


## CONTACT US

Please visit our website or call us with any inquiries:

Ph: 513-636-4685

Fax: 513-636-3861

[www.cchmc.org/DIL](http://www.cchmc.org/DIL)



## New Tests Now Available:

- Th17

## Tests Coming Soon:

- ALPS Biomarkers: soluble FAS-ligand (sFASL), IL-10, and IL-18

## Feedback:

We would like to hear from our Clients. We invite you to share your questions and comments with us. Feel free to send/fax/email your comments to us: Fax 513-636-3861; Email: [immunodeficiencies@cchmc.org](mailto:immunodeficiencies@cchmc.org)



## “Hat Day” in the DIL during the 2013 National Medical Laboratory Professionals Week.

(Clockwise left to right, Pat Adkins, Joyce Collett, Barb Wanstrath, Victor Lafay, Carrie Gifford and Kathryn Quinn.)

## NEW Fillable Requisition on the Website:

| GENERAL IMMUNOPHENOTYPING AND FUNCTIONAL TESTING   |   |  |  |
|--|---|--|--|
| <b>T Cell Immunophenotyping:</b><br><input type="checkbox"/> Lymphocyte Subsets<br><input type="checkbox"/> CD45 RA/RO<br><input type="checkbox"/> TCR V beta Repertoire<br><input type="checkbox"/> TCR alpha/beta/gamma/delta<br><input type="checkbox"/> Lymphocyte Activation Markers<br><input type="checkbox"/> iNKT cell enumeration  | <b>T cell Functional Studies :</b><br><input type="checkbox"/> Mitogen Stimulation<br><input type="checkbox"/> Antigen Stimulation<br><input type="checkbox"/> CD40L/ICOS<br><input type="checkbox"/> CTL Function<br><input type="checkbox"/> Cytokines, Intracellular Panel<br><input type="checkbox"/> pSTAT5  | <b>B Cell Immunophenotyping:</b><br><input type="checkbox"/> Lymphocyte Subsets<br><input type="checkbox"/> B cell panel<br><b>NK cell Functional Studies:</b><br><input type="checkbox"/> NK cell Function<br><input type="checkbox"/> CD107a (Note: cannot be accepted on Fridays) | <b>Cytokines &amp; Plasma Proteins</b><br><input type="checkbox"/> BAFF<br><input type="checkbox"/> Cytokine Panel: Plasma or CSF<br><input type="checkbox"/> Neopterin: Plasma or CSF<br><input type="checkbox"/> Soluble CD163<br><input type="checkbox"/> Soluble IL2 Receptor/s CD25 CD25<br>Cytokine Panel: IL-1B,2,4,5,6,8,10<br>IGN-g,TNF-a, and GM-CSF |
| DISEASE TARGETED TESTS   |   |  |  |
| <b>Disorders of Lymphoproliferation and Autoimmunity:</b><br><i>Autoimmune Lymphoproliferative Syndrome (ALPS)</i><br><input type="checkbox"/> ALPS Panel<br><input type="checkbox"/> Apoptosis (Fas-Mediated)<br>Note: Sample must be <24 hrs old and received Thursday AM<br><i>Immune dysregulation, Polyendocrinopathy, Enteropathy, X-linked Syndrome (IPEX)</i><br><input type="checkbox"/> Foxp3<br><br><b>Hemophagocytic Disorders:</b><br><i>Hemophagocytic Lymphohistiocytosis (HLH)</i><br><i>X-linked Lymphoproliferative Diseases (XLP)</i><br><input type="checkbox"/> NK function<br><input type="checkbox"/> CTL function<br><input type="checkbox"/> Perforin/Granzyme B<br><input type="checkbox"/> Soluble IL2 Receptor/sCD25<br><input type="checkbox"/> CD107a (Note: cannot be accepted on Fridays)<br><input type="checkbox"/> Soluble CD163<br><input type="checkbox"/> SAP/XLP1 (if male patient)<br><input type="checkbox"/> XIAP/XLP2 (if male patient)<br><input type="checkbox"/> Neopterin<br><input type="checkbox"/> iNKT<br><br><b>Immune Reconstitution Studies:</b><br><input type="checkbox"/> Lymphocyte Subsets<br><input type="checkbox"/> Mitogen Stimulation<br><input type="checkbox"/> NK Function<br><input type="checkbox"/> B Cell Panel<br><input type="checkbox"/> CD45 RARO<br><input type="checkbox"/> Foxp3 (T-reg enumeration) | <b>Combined Immune Deficiencies:</b><br><i>X-linked SCID</i> <input type="checkbox"/> CD127/132 <input type="checkbox"/> pSTAT5<br><i>Jak3 deficient SCID</i> <input type="checkbox"/> pSTAT5<br><i>IL-7R deficient SCID</i> <input type="checkbox"/> CD127/132<br><i>MHC Class I or II deficiency</i> <input type="checkbox"/> MHC Class I & II<br><i>Zap70 deficient SCID</i> <input type="checkbox"/> Zap70<br><br><b>Primarily Humoral Disorders:</b><br><i>Hyper IgM Syndrome (HIGM)</i><br><input type="checkbox"/> CD40L/ICOS<br><i>Combined Variable Immune Deficiency (CVID)</i><br><input type="checkbox"/> B Cell Panel<br><input type="checkbox"/> CD40L/ICOS<br><input type="checkbox"/> BAFF<br><br><b>Other Well-Defined Disorders:</b><br><i>Wiskott-Aldrich Syndrome (WAS)</i><br><input type="checkbox"/> WASP<br><input type="checkbox"/> WASP Post-Transplant Monitor Panel<br><i>Paroxysmal Nocturnal Hemoglobinuria (PNH)</i><br><input type="checkbox"/> PNH Screen (CD59) | <b>Neutrophil Disorders:</b><br><i>Chronic Granulomatous Disease (CGD)</i><br><input type="checkbox"/> Oxidative Burst (DHR)<br><i>Leukocyte Adhesion Deficiency (LAD)</i><br><input type="checkbox"/> Adhesion Markers (CD18/CD11b)   | <b>Other Testing and Services:</b><br><input type="checkbox"/> Lymphocyte Activation Markers (Basiliximab monitor)<br><input type="checkbox"/> CD52 (Alemtuzumab monitor)<br><input type="checkbox"/> CD64 (upregulation on neutrophils)<br><input type="checkbox"/> EBV Immortalized Line<br><input type="checkbox"/> Other:                                  |
| <b>Hyper IgE Syndrome:</b><br><input type="checkbox"/> Th17 Enumeration  |   |  |  |

### IN THE UPCOMING ISSUE:

- ALPS Biomarkers

