

**Heart Institute Diagnostic Lab**

CAP#: 1667801

CLIA#: 36D0656333

Phone: (513) 803-1751

Fax: (513) 803-1748

Email: HeartDx@cchmc.org

## Familial Hypercholesterolemia

### *LDLR Full Gene Sequencing* *APOB (R3500W and R3500Q) Sequencing*

Familial Hypercholesterolemia (FH) is characterized by elevated LDL cholesterol levels which can lead to atherosclerotic plaque depositions in the arteries. FH is associated with a markedly increased risk of coronary artery disease at a young age. Deposits of cholesterol can also be seen in tendons (xanthomas) or around the eyes (xanthelasmas). Individuals with FH commonly have coronary heart disease, which can be associated with angina, myocardial infarction, and stroke.

The prevalence of heterozygous FH is estimated to be 1:200-500. Untreated individuals have a 20-fold increased risk for coronary heart disease. The diagnosis of FH can be confirmed by the presence of a pathogenic variant in genes associated with FH. Mutations, including deletions and duplications, in *LDLR* can be identified in 60-80% of individuals with FH. Pathogenic variants in *APOB* account for about 1-5% of FH cases, with R3500W and R3500Q mutations being the most common. FH is an autosomal dominant condition; however, homozygous FH can occur with a more severe phenotype.

### Shipping Instructions

Please enclose a test requisition form with sample. All information must be complete before sample can be processed. Samples may be shipped at room temperature by overnight Federal Express to arrive Monday through Friday.

### Ship To:

Cincinnati Children's  
Hospital Medical Center  
Attn: Heart Institute Diagnostic Lab  
240 Albert Sabin Way,  
Room S4.381  
Cincinnati, OH 45229-3039

### Indication

Genetic testing for FH is indicated for individuals with clinical suspicion for FH.

Currently available tests include:

LDLR Full Gene Sequencing

LDLR Known Mutation Sequencing

APOB (R3500W/R3500Q) Analysis

## Methodology:

All 18 coding exons, as well as the exon/intron boundaries and portion of untranslated regions of the *LDLR* gene (NM\_000527.4) are amplified by PCR for *LDLR* full gene sequencing analysis. A region harboring the R3500 of the *APOB* gene is amplified by PCR for the *APOB* (R3500W and R3500Q) sequencing test. Genomic DNA sequences from both forward and reverse directions are obtained by automatic fluorescent detection using an *ABI PRISM® 3730 DNA Analyzer*. Sequence variants different from National Center for Biotechnology Information GeneBank reference are further evaluated for genetic significance. If a mutation is identified, a known familial mutation analysis will be available for additional family members.

## Sensitivity & Accuracy:

Greater than 99% of the sequencing mutations in exons 1-18 of *LDLR* and *APOB* R3500W and R3500Q mutations are detectable by this method. Sequencing does not detect large DNA rearrangements, deletions/duplications or low level mosaicism.

## References:

Goldberg AC, et al. Familial hypercholesterolemia: Screening, diagnosis and management of pediatric and adult patients: Clinical guidance from the national lipid association expert panel on familial hypercholesterolemia. *Journal of Clinical Lipidology*. 2011;5:S1-8.

*Identification and management of familial hypercholesterolaemia (FH)*. London; 2008.

Youngblom E, Knowles JW. Familial hypercholesterolemia. In: Pagon RA, Adam MP, Ardinger HH, Wallace SE, Amemiya A, Bean LJH, Bird TD, Dolan CR, Fong CT, Smith RJH, Stephens K, eds. *GeneReviews(r)*. Seattle (WA); 1993.

## Specimen:

Peripheral blood in EDTA tube

Adult: 3-5mL

Child: 3-5mL

Infant: 1-3mL

For other specimen types, please contact Amy Shikany at 513-803-3317

## Turnaround Time:

LDLR Full Gene Analysis 2-4 weeks

LDLR Known Mutation Analysis 1-2 weeks

APOB (R3500W/R3500Q) Analysis 1-2 weeks

## CPT Codes:

LDLR Full Gene Sequencing: 81406

LDLR Known Sequencing: 81403

APOB (R3500W/R3500Q) Analysis: 81479