

# Organotypic Raft Cultures

1

## Materials:

1. Cotton Pads - #740-E, 8 x 10 in., 25/pack. Cut into 1 x 1 in squares. Autoclave in jar.
2. Sterile Forceps - Autoclave in paper bag.
3. Transwell Inserts - Order from Fisher (Costar #3450)
4. Organogenesis Trays (BD Biosciences #355467)
5. Rat Tail Collagen Type 1 - Upstate Biotechnology Inc. #08-115.
6. C8:0 - 1,2-dioctanoyl-*sn*-glycerol (Sigma #D1912)

## Medias:

1. Fibroblast Media:  
F12 media  
10% FBS  
1X pen/strep
2. Keratinocyte Plating Media:  
F media (3:1 F12/DMEM) = 0.66 mM Ca<sup>2+</sup>  
0.5% FBS  
1X (hydrocortisone, cholera toxin, insulin, adenine, pen/strep)  
1.88 mM Ca<sup>2+</sup> (add 305  $\mu$ l 1M CaCl<sub>2</sub> to 250 ml of media)
3. Cornification Media 1:  
F media (3:1 F12/DMEM) = 0.66 mM Ca<sup>2+</sup>  
5% FBS  
1X (hydrocortisone, cholera toxin, insulin, adenine, pen/strep)  
1.88 mM Ca<sup>2+</sup> (add 305  $\mu$ l 1M CaCl<sub>2</sub> to 250 ml of media)
4. Cornification Media 2:  
F media (3:1 F12/DMEM) = 0.66 mM Ca<sup>2+</sup>  
2.5% FBS  
1X (hydrocortisone, cholera toxin, insulin, adenine, pen/strep)  
1.88 mM Ca<sup>2+</sup> (add 305  $\mu$ l 1M CaCl<sub>2</sub> to 250 ml of media)
5. 10X F12 media:  
Dissolve 1 packet Ca<sup>2+</sup> free F12 in 90 mls ddH<sub>2</sub>O  
Add 1.176 g of sodium bicarbonate (NaHCO<sub>3</sub>)  
pH to 7.2  
Add 23  $\mu$ l of 1M CaCl<sub>2</sub>  
Filter sterilize and store at -20°C

## Preparation of Collagen Raft:

2

1. Trypsinize fibroblasts. Count and make final concentration  $7.5 \times 10^5$  cells/ml of Fibroblast media.
2. Prepare collagen premix (for 25 ml):
  - a. 10X F12 media 2.5 ml
  - b. 10N NaOH 6  $\mu$ l
  - c. Pen/strep 250  $\mu$ l
  - d. FBS 2.5 ml
  - e. collagen solution 20 mlMix carefully. Avoid bubbles.
3. Plate 1 ml of collagen premix per Transwell insert. Tilt to coat bottom. Allow to gel ~5 minutes in the hood.
4. Add 600  $\mu$ l of  $7.5 \times 10^5$  cells/ml of fibroblasts ( $4.5 \times 10^5$  cells) to remaining collagen mixture. Mix thoroughly. Avoid bubbles.
5. Layer 2.6 ml of fibroblast/collagen solution over each 1 ml gel. Allow to gel ~30 minutes in incubator.
6. Add ~20 mls Fibroblast media to outer well (or enough to submerge the collagen raft. Note: Once raft starts to float, you have enough media. The media will rise through the transwell membrane and through the raft itself.). Incubate in 5% CO<sub>2</sub>, 37°C incubator.
7. Gels should contract to appropriate shape within 4-7 days. You do not need to change the media during this time.

## Plating of Keratinocytes:

3

### Day 0:

1. Aspirate media from outer well. Carefully remove media from inner Transwell by using a pipette. Layer 150  $\mu$ l of  $1.4 \times 10^6$  cells/ml ( $7 \times 10^5$  cells) onto collagen raft.
2. Incubate for 2 hours to allow cells to attach. Carefully add ~20 ml Keratinocyte Plating media to outer well.

### Day 2:

1. Aspirate media from outer well. Carefully remove media from inner Transwell by using a pipette.
2. Feed cells with ~20 mls Keratinocyte Plating media to outer well.

### Day 4:

1. Aspirate media from outer well. Carefully remove media from inner Transwell by using a pipette.
2. Sterily lift transwell insert and place 2 cotton pads on the risers of each outer well of the 6 well Organogenesis plate. Place transwell insert on the cotton pads.
3. Add 10-11 mls Cornification Media 1 containing 10  $\mu$ M C8:0 (to induce productive HPV life cycle) to each outer well.

Feed cells every other day with 10 mls Cornification Media 1 containing 10 $\mu$ M C8:0. At about Day 9, appearance of skin equivalent should be dry. Full stratification occurs by Day 15. Full HPV life cycle peaks ~Day 16. For long experiments, switch to Cornification Media 2 containing C8:0 for feeding post-Day 15.

## BrdU Labeling:

1. Aspirate media from outer well. Carefully remove media from inner Transwell by using a pipette.
2. Add 10 ml Cornification Media 1 containing 10  $\mu$ M C8:0 and 10  $\mu$ M BrdU to each outer well.
3. After 8 hours, fix in 4% formalin (see Elsa's Raft Fixation/Preparation for Histology Protocol).